

Products by a Dentist, for Dentists

ORDER FROM YOUR DENTAL SUPPLY COMPANY



NEW THE VACU-MIXER AUTOMATED VPS MIXING MACHINE STARTER KIT

- Vacu-Mixer and Charger / Adapter
- 50 Dynamic Mixing Tips
- 50 Harmony Posterior Trays
- 10 Vaccu-sil 120ml Cartridges
- 50 Mojo Complete Syringes
- 4 VacuCord Gingival Retraction Cord
- Choice of Regular set or Fast Set



Heavy body Regular Set Heavy body Fast Set Light body Regular Set Light body Fast Set

#1 RATED VACCU-SIL SUPER HYDROPHILIC VPS IMPRESSION MATERIAL

NEW THE VACU-MIXER AUTOMATED VPS MIXING MACHINE



NEW VACCU-SIL 120 ML CARTRIDGES
Heavy Regular | Heavy Fast

NEW DYNAMIC MIXING TIP



NEW VacuCord Gingival Retraction Cord



NEW VACCU-SIL BITE REGISTRATION
50 ML CARTRIDGES



MYLAR MATRIX ROLL
9 Microns thin



THE SAND TRAP



HO SECTIONAL MATRIX
ORIGINAL, CONTOURED, AND
CONTOURED DEEP/WIDE



THE MOJO SYRINGE
Excellent storage and delivery system for light body or mono-phase VPS.
Ready to load.



NEW THE MOJO COMPLETE
The Mojo syringe pre-filled with Vaccu-sil light body VPS.
Regular or Fast set



HARMONY ANTERIOR
DUAL ARCH TRAY
WITH FACIAL PLANE RELATOR FEATURE



FACIAL PLANE RELATOR



HARMONY POSTERIOR
DUAL ARCH TRAY

HO DENTAL COMPANY FAMILY OF PRODUCTS

ORDER FROM YOUR DENTAL SUPPLY COMPANY

The Harmony Anterior Dual Arch Tray

The Harmony Anterior Dual Arch Tray will capture the prepared teeth, the opposing arch, the bite registration, and the patient's facial midline, in one simple step. The Harmony Anterior Dual Arch Tray can be used without alignment sticks for single units, or with them for multiple units. No more slanted crowns or veneers.

Harmony Anterior Tray HTT103 40 Trays / 40 sticks /box

The Harmony Posterior Dual Arch Tray

The narrow stainless steel bar at the back is only 2.6 mm in diameter, preventing tuberosity impingement, while still giving it rigidity equal or better than any other disposable posterior dual arch trays.

Harmony Posterior Tray HTT104 50/box

Vaccu-sil VPS Impression Material

In an extensive study done in 2014 by an Utah-based independent dental products evaluation organization, Vaccu-sil was rated number one in the world amongst all VPS impression materials based on value and performance.

Vaccu-sil Heavy Regular	Vac-50HR	4 50ml cartridges
Vaccu-sil Heavy Fast	Vac-50HF	4 50ml cartridges
Vaccu-sil Light Regular	Vac-LR-Cart	4 50ml cartridges
Vaccu-sil Light Fast	Vac-LF-Cart	4 50ml cartridges

NEW Vaccu-sil Bite Registration

Like our highly rated Vaccu-sil heavy and light body VPS materials, our Vacu-Bite stands up to the major brands and yet is the lowest cost on the market.

Vaccu-sil Bite Registration Vacu-Bite 2x50ml cartridges

NEW VacuCord Gingival Retraction Cord

A knitted retraction cord used for rapid tissue displacement and absorption of fluids at the sulcus, in preparation for a dental impression. Sizes 00, 0, 1, and 2.

Individual	VacuCord 00	1 unit (100 inch)
	VacuCord 0	1 unit (100 inch)
	VacuCord 1	1 unit (100 inch)
	VacuCord 2	1 unit (100 inch)
All Four Sizes (00, 0, 1, 2)	VacuCord 4	4 units, (100 inches each)

NEW The Vacu-Mixer Automated VPS Mixing Machine

The Vacu-Mixer is hand held, portable, and rechargeable. Weighing less than two pounds, the Vacu-Mixer holds a 120 ml cartridge instead of the usual 50 ml. The larger volume cartridge lowers the cost of materials, has less waste, most importantly, it delivers a perfect mix each and every time.

Vacu-Mixer Machine Vacu-Mixer 1 Machine

NEW The Vacu-Mixer Automated VPS Mixing Machine Starter Kit

Excellent Value! Start enjoying the instant benefits of our hand-held automix unit along with over \$400 of free supplies! Choice of regular or fast set formulas.

Vacu-Mixer Machine	Vacu-Mixer	1 Machine
Adapter/Charger	Adapter/Charger	Adapter/Charger
Dynamic Mixing Tips	Vaccu-Mix Tip	50 pc / bag
Vaccu-sil 120ml Cartridges	Regular or Fast	10x120ml Cart/box
Harmony Posterior Tray	HTT104	50 / box
Mojo Complete	Regular or Fast	50 syringes/box
VacuCord	VacuCord 4	Sizes 00,0,1,2 / 100 In.

For Regular Set order part #: Vacu-Mix-Kit-R

For Fast Set order part #: Vacu-Mix-Kit-F

NEW Dynamic Mixing Tips

The Dynamic Mixing Tips are compact and save material usually wasted in larger mixing tips.

Dynamic Mixing Tip Vacu-Mix Tip 50pc/bag

NEW Vaccu-sil Heavy Body VPS Impression Material 120ml Cartridges

Vaccu-sil Heavy Reg	Vac-120HB-R	10X120ml Cart/box
Vaccu-sil Heavy Fast	Vac-120HB-F	10X120ml Cart/box

The Facial Plane Relator

With Ho Dental's Facial Plane Relator, you can easily record an accurate bite registration and the exact position of the desired midline in relation to a horizontal plane in just 60 seconds. The Facial Plane Relator will accurately communicate the patient's vertical and horizontal midlines to the laboratory technician.

Facial Plane Relator (per unit)	FPR102	100/box
Facial Plane Relator (per dozen)	FPR102/12	12/box
Facial Plane Relator (24 pack)	FPR102/24	24/box

The Mojo Syringe (ready to fill)

The material wasted inside a standard mixing tip (up to 2.4 CC) is enough to do 3-4 more crowns with our syringe. The extra-narrow tip enables you to apply the impression material right at the sulcus, circumferentially, without having to lift the tip, thus, minimizing air bubbles and retakes. The ergonomic design allows easy access even to the distals of second molars.

Mojo Syringe Mojo101 100/Box

NEW The Mojo Complete (replaces the BFC Complete)

The Mojo Syringe prefilled Vaccu-sil light body VPS impression material, for the ultimate in convenience. Available in regular or fast setting formulas. Eliminate cross-contamination and messy clean ups. Replaces the BFC Complete with more storage and improved controls.

Light Body Regular Set	Mojo-CR	50 syringes/box
Light Body Fast Set	Mojo-CF	50 syringes/box

Ho Sectional Matrix

From the creator of the HO Band, the Ho Sectional Matrix is so perfectly sized and contoured that it's all you need to complete your restorations. And, like the original Ho band, the Ho Sectional Matrix is 0.001 inch thin. Available in deep/wide (D/W) and contoured.

HO Section Original Matrix	HSM201	100/box
HO Sectional Contoured D/W Matrix	HSM204	100/box
HO Sectional Contoured Matrix	HSM206	100/box

Mylar Matrix Roll

Our Mylar is ultra-thin. At a thickness of just nine microns, it can be used during cementation and removed afterwards without causing a gap. Uses include cementation of porcelain veneers and other ceramic restorations.

Mylar Matrix Roll MLR202 300ft roll

The Sand Trap

The Sand Trap is a patented clear plastic device designed to confine and suction away aluminum oxide particulates during intra-oral or chair side sandblasting.

Sand trap (bag of 100 units)	HST205/100	100/bag
Sand trap (bag of 12 units)	HST205/12	12/bag